

[Protocol 4.01] Testing for remaining Sendai virus_CytoTune2.0

Author: Jeong-Eun Lee
No tags associated

Created: 08.12.2021 11:01
Last modified: 12.12.2023 09:48

No custom dates added

Protocol 4.01 Testing for remaining Sendai virus_CytoTune2.0
Version: 1.0 (12.12.2023)

Media and Reagents:

RNeasy Plus Mini Kit (Qiagen, 74134)
SuperScript III First-Strand Synthesis System
(Thermo Fisher,18080051)
Taq polymerase (Thermo Fisher, 10342020)
100 mMdNTP Set (Thermo Fisher,10297-018)
Agarose (Thermo Fisher,15510-027)
GelRed (Biotium, 41003)
10x BlueJuice Gel Loading Buffer (Thermo Fisher,10816015)
100 bp Ladder (Thermo Fisher,15628-019)

Materials and Equipment:

0,2 ml PCR tubes
1,5 ml tubes
PCR Cyclor
Electrophoresis chamber
Nanodrop photometer

1. Introduction and Purpose

This protocol describes a RT-PCR assay to test for remaining Sendai virus genome and transgenes in iPSC cultures after reprogramming with CytoTune-iPS 2.0 Reprogramming Kit (Thermo Fisher, A16517).

2. RNA isolation and reverse transcriptase

- Isolate RNA using the RNeasy Plus Mini Kit following manufactures instructions. Isolate RNA from two separate wells of a six well plate using separate columns for each wells.
- Determine the amount and quality of RNA using a Nanodrop photometer.
- Use 1 µg RNA with SuperScript III First-Strand Synthesis System for reverse transcription following the manufacturer's protocol (The half amount of the recommended protocol is sufficient).
- Dilute cDNA after reverse transcription 1:10 (20 µl + 180 µl) with water.

3. Procedure Sendai Virus PCR

- The primer pairs for CytoTune Sendai virus detection and for a control housekeeping gene (Hu18SRNA) are listed in the table below.

Target	Internal primer name	Primer sequence	T _m °C	Product size [bp]
SeV	SeV_F	GGA TCA CTA GGT GAT ATC GAG C*	55	181
	SeV_R	ACC AGA CAA GAG TTT AAG AGA TAT GTA TC*	55	
Klf4	SeV_Klf4_F	TTC CTG CAT GCC AGA GGA GCC C	55	410
	SeV_Klf4_R	AAT GTA TCG AAG GTG CTC AA*	55	
cMyc	SeV_cMyc_F	TAA CTG ACT AGC AGG CTT GTC G*	55	532
	SeV_cMyc_R	TCC ACA TAC AGT CCT GGA TGA TGA TG	55	
KOS	SeV_KOS_F	ATG CAC CGC TAC GAG TGA GCG C	55	528
	SeV_KOS_R	ACC TTG ACA ATC CTG ATG TGG	55	
Hu18SRNA	Hu18SRNA_F	GTA ACC CGT TGA ACC CCA TT	50	151
	Hu18SRNA_R	CCA TCC AAT CGG TAG TAG CG	50	

*Primer contains SeV genome sequences. Pairing of these primers with transgene-specific primers allows specific detection of transgenes carried by the CytoTune Sendai reprogramming vectors.

- Use 2 µL of diluted cDNA for Sendai Virus PCR.
- PCR is set up as follows:

PCR ingredient	µl per reaction
10x PCR buffer	2
dNTP-Mix (10 mM each)	0.8
50mM MgCl2	0.6
Forward primer (100 pmol /µl)	0.5
Reverse primer(100 pmol /µl)	0.5
Water	13.5
Taq polymerase (5 U/µl)	0.1
diluted cDNA	2
total	20 µl

Note:
Negative control 2 µl of diluted cDNA from a control PSC line not derived using Sendai virus
Positive control 2 µl of diluted cDNA from generated iPS (passage 2 or 3) after reprogramming

4. Close and spin all PCR tubes briefly, load the PCR cyclcr and start the program.

Cycler program:

1 cycle	94 °C	3 min
35 cycles	95 °C	30 sec
	55 °C	30 sec
	72 °C	30 sec
1 cycle	72 °C	10 min
	4 °C	∞

5. Cast a 2 % standard agarose gel with DNA stain GelRed 5µl/100 ml.
6. Pipette 2 µl of 10x BlueJuice to each PCR reaction. Use 10 µl of PCR reaction for gel.
7. Use 1 µl of 100 bp Ladder mixed with 1 µl of 10x BlueJuice and 8 µl of water as size and running control.
8. Run agarose gel at 80 -100 V for 20 - 40 minutes.

Relevant applicable documents:
Manual CytoTune-iPS 2.0 Reprogramming Kit (Thermo Fisher, A16517)
Manual SuperScript III First-Strand Synthesis System (Thermo Fisher,18080051)
Manual RNeasy Plus Mini Kit (Qiagen, 74134)