

[Protocol 1.03.1] Reconstitution, aliquoting and use of Y26732 ROCK inhibitor

Author: Jeong-Eun Lee
No tags associated

Created: 20.12.2021 18:56
Last modified: 11.08.2022 14:19

No custom dates added

Protocol 1.03.1 Reconstitution, aliquoting and use of Y26732 ROCK inhibitor
Version: 1.0 (20.12.2021)

Media and Reagents:

ROCK inhibitor (Y-27632; Stem Cell Technologies, 72305)
Distilled sterile water

Alternatively in DMSO, ready to use:

ROCK inhibitor (Y-27632 2HCl; Selleck Chemicals SEL-S1049-10mM pre-solved)

Materials and Equipment:

Sterile syringe filter (low protein-binding!)
0,5 ml Eppis

1. Introduction and Purpose

This protocol describes the storage of Y-27632 ROCK inhibitor and the preparation of media containing Y-27632 ROCK inhibitor. Y-27632 is a $\geq 95\%$ pure synthetic compound with a molecular weight of about 325 g/mol. It is light-sensitive and hygroscopic.

2. Reconstitution and aliquoting

Centrifuge the vial (Stem Cell Technologies, 72305) prior to opening. Reconstitute in sterile water.

Reconstitute Y-27632	in sterile water	Concentration of solution
10 mg	5,912 ml	5 mM

Filter sterile using a low protein binding filter and prepare aliquots of 200 μ l. Store these at -20°C.

Alternatively, 10 mM ROCKi in DMSO, ready to use (Selleck Chemicals):

Prepare aliquots of 50 -100 μ l and store these at -20°C or at -80°C for long-term storage up to two years.

The location of aliquots: **TBA**

3. Preparation of medium containing ROCK inhibitor

Thawed aliquots of Y-27632 ROCK inhibitor can be stored at 4°C for up to two weeks.

The final concentration of Y-27632 ROCK inhibitor in the culture medium should be 10 μ M. Therefore, the stocks will be used 500x (5mM) or 1000x (10mM) diluted.

Following table shows the apposite volume of Y-27632 ROCK inhibitor for different volumes of medium.

Volume Medium [ml]	Volume ROCKi 5mM (Stem Cell Tech) [μ l]	Volume ROCKi 10mM (Selleck Chem) [μ l]
5	10	5
12	24	12
50	100	50
500	1000	500

Relevant applicable documents:

Protocol 1.03.0 Passaging of iPSC into single cells using TrypLE/Accutase
Protocol 1.04.0 Sorting of differentiated cells using MACS