

[Protocol 1.01.1] Freezing of iPSCs using Bambanker freezing medium

Author: Jeong-Eun Lee
No tags associated

Created: 11.08.2022 12:06
Last modified: 07.09.2022 14:28

No custom dates added

Protocol 1.01.1 Freezing of iPSCs using Bambanker freezing medium
Version: 1.0 (11.08.2022)

Media and Reagents:

Bambanker™ (Nippon Genetics, BB01)
0,5 M EDTA (Thermo Fisher, 15575-020)

Materials and Equipment:

Aspiration Pump
Cell scraper
Freezing container (Mr. Frosty) filled with 100% 2-propanol (pre-chilled)
Pre-labeled Cryo vials (NUNC, 479-6843)

1. Introduction and Purpose

This protocol describes the methodology for cryopreservation of induced pluripotent stem cells (iPSCs) using Bambanker™ freezing medium.

2. Procedure Cryopreservation

- Cells must be cryopreserved when they are in their log phase of growth (60 - 70 % confluence).
- Use a microscope to visually identify regions of differentiation and remove regions of differentiation by scraping with a pipette tip. This selection should not exceed 20% of the well if the culture is of high quality.
- Aspirate medium from the culture and rinse once with appropriate volume of PBS (see table below).

Culture Vessel	Surface Area [cm2]*	Volume 0.5 mM EDTA [ml]	Volume Bambanker TM[ml]
24 well plate	2	0.25	0.25
12 well plate	4	0.5	0.5
6 well plate	10	1	1
T25 flask	25	3	3
T75 flask	75	6	6

*Exact sizes are manufacturer/type dependent

- To lift cells from the tissue culture vessel add appropriate volume of 0.5 mM EDTA to the culture (see table above) and incubate for 3 - 5 minutes at RT.

Note: Alternatively the incubation can occur at 37°C on heating plate in the IVF bench. In this case the incubation time will be reduced to 2 - 3 minutes.

- Check under the stereo microscope, if cell colonies start to show holes throughout. Care should be taken not to incubate until the colonies detach from the plate.
- Aspirate EDTA gently and add appropriate volume of Bambanker™ freezing medium (see table above) to the cells.
- Gently tap the plate 5 - 10 seconds to dislodge the cells from the plastic and use a 5 ml pipette to fill the cell suspension (ideally 1 ml) into labeled cryo vials. Do not aspirate more than 2 times to avoid breaking the cell clumps into single cells.
- Note:** Alternatively the cells can also be scraped off with a cell scraper. Take care to keep the clumps as big as possible!
- Place the cryo vials immediately into the freezing container and place at -80°C overnight.
- Transfer the cells to a liquid N2 tank the next day.

Appendix

General points to consider:

- Wells are 60 - 70% confluent at time of cryopreservation.
- Cryopreservation should be done approximately 1 day before the cells are ready to passage.
- Before cryopreservation, pluripotent stem cells should be of high quality (primarily undifferentiated with less than 20% of the cells being differentiated).
- Pluripotent stem cells will have improved survival following thawing if cryopreserved as large clumps.
- 1 ml of cryopreservation media should be used for every well of a 6-well plate being frozen. However, if the wells are at low density (less than 50% confluent), 1 ml of cryopreservation media may be used for every 2 wells.

Label of cryo vials as follows and place on ice:

- Cell type/cell line name
- Passage
- Freeze date
- Your Initials

Relevant applicable documents:

Protocol 1.01.2 Thawing of hPSC from Bambanker™ freezing medium

