

[Protocol 2.02] NGN2 protocol

Author: Jeong-Eun Lee
No tags associated

Created: 25.05.2022 19:00
Last modified: 22.02.2023 10:30

No custom dates added

Protocol 2.02 NGN2 protocol
Version: 1.0 (25.05.2022)
According to the [Südhof paper](#)

Media and Reagents:

DMEM/F-12, HEPES (11330-032, Thermo Fisher)
Neurobasal Medium (21103-049, Thermo Fisher)
N2 supplement (100x) (17502-048, Thermo Fisher)
MEM Non-Essential Amino Acids Solution (100X) (11140-035, Thermo Fisher)
B27 supplement (50x) (17504-044, Thermo Fisher)
GlutaMAX Supplement (35050061, Thermo Fisher)
hBDNF (450-02, peprotech)
hNT-3 (267-N3-025, R&D systems)
Laminin (L2020, Sigma-Aldrich)
doxycycline (D3447-500MG, Sigma-Aldrich)
puromycin (4089, tocris)
Virus
DPBS, no calcium, no magnesium (14190250, Thermo Fisher)
DPBS, calcium, magnesium (14040091, Thermo Fisher)

Materials and Equipment:

6/12/24 well plates
10ul/100ul/200ul/100ul pipettes & pipette tips
5ml/10ml pipettes
Cell counting equipment
15/50 ml Falcon tube
1,5 ml sterile Eppis

1. Introduction and Purpose

This protocol describes the protocol converting human iPSCs into functional iN cells with nearly 100% yield and purity in less than 2 weeks by forced expression of a single transcription factor, NGN2. The resulting ES-iN or iPS-iN cells exhibit quantitatively reproducible properties independent of the cell line of origin, form mature pre- and postsynaptic specializations.

2. Protocol

Day -2: seed cells for neuronal differentiation

- Aspirate medium
- Add 1ml (6er well) /3ml (10 cm) Accutase
- Incubate 6-8 min at 37°C
- Resuspend and add 10 ml E8 home
- Spin: 1000 rpm for 3-5 min
- Resuspend in E8 home + ROCK inhibitor
- Count
- Dilute cell suspension
- Aspirate matrigel
- Seed cells, e.g. 6-well: seed 200.000 cells/well (next day around **40 %** confluent)

Day -1: virus infection

- Mix E8 and viruses
- Viruses: 4ul NgN+GFP + 3ul rtTA per 6-well (Define the optimal virus titer if necessary)

Day 0: medium change and dox activation

- Aspirate medium and replace with F12-N2 medium + **3 µg/ml of dox**
- Add dox always fresh to the media! (The dox aliquot)
- The thawed dox aliquot can be stored at 4°C for one week

Day 1: medium change

- Aspirate medium and replace with F12-N2 medium + **3 µg/ml of dox**

Day 2: puromycin selection

- Aspirate medium and replace with F12-N2 medium + 3 µg/ml of dox
 - add 0.8 ug/ml puromycin (stock: 10 mg/ml 1:12.500) to F12-N2 medium
- Day 3: adding the mouse glia
- Detach glia using trypsin, spin, resuspend in NB-B27, 1x T75 for 72 x 24 wells / 14x 6wells-suck up medium, add medium containing glia enhancing synapse formation
- Day 5 - 14: media change
- change 50% medium NB-B27 every other day

3. Media preparation

F12-N2 media:

Component	Stock concentration	Final concentration	50 ml
DMEM/F12			48.9 ml
N2 supplement	100 x	1 x	0.5 ml
NEAA	100 x	1 x	0.5 ml
hBDNF	10 µg/ml	10 ng/ml	50 µl
hNT-3	50 µg/ml	10 ng/ml	10 µl
laminin	1 mg/ml	0.2 µg/ml	10 µl
doxycycline	2 mg/ml	3 µg/ml	add directly before use

NB-B27 media:

Component	Stock concentration	Final concentration	50 ml
Neurobasal medium			48.5 ml
B-27 supplement	50 x	1 x	1 ml
GlutaMAX	100 x	1 x	0.5 ml
hBDNF	10 µg/ml	10 ng/ml	50 µl
hNT-3	50 µg/ml	10 ng/ml	10 µl
laminin	1 mg/ml	0.2 µg/ml	10 µl
AraC	4mM	2µM	Add only at day 4 (1:2000)
doxycycline	2 mg/ml	3 µg/ml	add directly before use

4. Confirmation of neuronal differentiation
- Immunofluorescence imaging - MAP2 & Tuj1 (pan neuronal markers)
 - qPCR